

Surgical Led Light (Portable)

SMart-M35

Tips

- 1. Double power supply system and built-in battery will start work within 0.2s when city power is off.
- 2. Intelligent charging mode will monitor battery and show its working state.
- 3. Average LED life ≥60000 hours
- 4. Control keypad: light held adjustment; ON/OFF; color changes; brightness
- 5. One button to change from ENDO, microwave and R9 mode
- 6. Spot diameter ≥300mm
- 7. Imported Germany gas spring for durable and exact positioning.



SPECIFICATIONS:

Standard Configuration				
	No.	Item	Brand	Quantity
	1	Major Dome	SHINVA	1 pc
	2	Main Body	SHINVA	1 pc
	3	Power Cord	SHINVA	1 set
	4	Light Source	Germany OSRAM	1 set
	5	Power Switch	MEANWELL	2 pcs
	6	Spring Arm	Germany ONDAL	2 pcs
	7	Rotator	Germany ONDAL	1 set
	8	Dust Cover	SHINVA	1 set
	9	Sterilized Hand Grip	SHINVA	4 sets
Technical Configuration		100	7.	
Center Illumination:	80000≤Ec≤160000			
Color Temperature Tc:	3500 ~4500K			
Illumination beam depth:	≥1200mm			
Deep Cavity Illumination Rate:	≥90%			
Color Rendering Index Ra:	96≤Ra≤100			
Irradiance W/m²:	≤500			
Irradiation Ratio mW/ (m2•lx):	≤4			
Built-in Camera (Optional)				
Brand:	Haikang			
Pixel:	2380000			
Optical Zoom:	10 times			
Aperture:	F=1.8 to 3.4			
Censor Type:	CMOS			
Censor Dimension:	1/2.8 inch			
Focal Length:	10mm (wide end), 1000mm (far end)			
Control Way:	Camera can be operated from surgical light control system, and there is no need to prepare separate camera control box.			
Camera System:	Meet shooting needs			

Tel: 8373-0097; 8899-3433 Telefax: 8374-3997 Email: sales.bajpharma@gmail.com

